



Life Sciences

Better Lives. Better Planet.SM



fortéBIO[®]
A Division of Pall Life Sciences

Revolution in Interactions Analysis

Attila Aranyos Ph.D.
2016 May

This presentation is the copyright work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer

© 2015 Pall Corporation



Life Sciences

fortéBIO[®]
A Division of Pall Life Sciences

Better Lives. Better Planet.SM



Bio Layer Interferometry Label Free Revolution in Protein Quantitation and Real- time Kinetics

- Working Principles
- The BLI instrumentation platform
- Difference / Complementarity with SPR
- Key applications with scientific examples


Attila Aranyos, Ph.D

This presentation is the copyright work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer

© 2015 Pall Corporation



Label-free , Real-time Techniques




Mainstream

- BLI (Bio Layer Interferometry)
- SPR (Surface Plasmon Resonance)

Others


- QCM (Quartz Crystal Microbalance)
- SAW (Surface Acoustique Wave)



- They differ in detection method, but they all provide « Sensograms »
- They largely differ in speed, throughput, sample compatibility, usability



FortéBio - a Division of Pall Life Sciences*



1990: Biacore

2006 : GE


- Historical leader
- Global Market share leader

2003: Other SPR Systems

- Never as good as Biacore

*Pall corporation is part of Danaher Corporation since Aug 2015

FortéBio - a Division of Pall Life Sciences



- What do we Propose:** Proprietary bio-layer interferometry (BLI) technology for label-free bio sensing (Immunoassay, kinetics and quantitation) as a credible alternative to SPR
- Key to success:** SPR-quality data without the inconveniences the SPR systems : no fluidics, more channels more flexibility, simpler to use
- Success factors 2015:**
 - 9th consecutive year of 25%+ growth
 - 1200 users in user meetings
 - Over 1500 systems in total , 15 million biosensors
 - 300+ publications, in respected papers including Nature
 - Multi-system industrial accounts : 2 out 3 users

5

PALL Life Sciences



Since the first Octet system shipment, Nearly 1200 articles published in peer reviewed journals

jbc THE JOURNAL OF BIOLOGICAL CHEMISTRY

CVI CLINICAL AND VACCINE IMMUNOLOGY
Published Monthly by the American Society for Microbiology
2006-PRESENT

nature structural & molecular biology

Science
AAAS

Protein Expression and Purification
ACADEMIC PRESS

PROTEIN SCIENCE

Cancer Research

VIROLOGY JOURNAL

JIM

PNAS

Cell Host & Microbe

BIOCHEMISTRY
including biophysical chemistry & molecular biology

PALL Life Sciences

1170
1000
800
600
400
200
0
Per year
Cumulative

Year	Per year	Cumulative
2005	7	7
2008	41	67
2010	101	225
2012	221	667
2014	309	1170

JMB
JOURNAL OF MOLECULAR BIOLOGY

Structure

ANALYTICAL BIOCHEMISTRY
Methods in the Biological Sciences
ACADEMIC PRESS

nature immunology



Since the first Octet system shipment, Nearly 1200 articles published in peer reviewed journals

jbc THE JOURNAL OF BIOLOGICAL CHEMISTRY

CVI CLINICAL AND VACCINE IMMUNOLOGY
Published Monthly by the American Society for Microbiology
2006-PRESENT

nature structural & molecular biology

Sci LETTER
BIOLOGY

doi:10.1038/nature13614

Protein Expression and Purification
ACADEMIC PRESS

LETTER
BIOLOGY

HSP70 sequestration by free α -globin promotes ineffective erythropoiesis in β -thalassaemia

Jean-Benoit Arlet^{1,2,3,4*}, Jean-Antoine Ribeil^{1,3,4,5*}, Flavia Guillem^{1,3,4}, Olivier Negre⁶, Adonis Hazoume^{7,8}, Guillaume Marcion^{7,8}, Yves Beuzard⁹, Michael Dussiot^{1,3,4}, Ivan Cruz Moura^{1,3,4,9,10}, Samuel Demarest¹¹, Israe Chauvet de Beauchene^{11,12}, Zakia Belaid-Choucair^{1,3,4}, Margaux Sevin^{7,8}, Thiago Trovati Maciel^{1,3,4,9,10}, Christian Auclair^{11,12}, Philippe Leboulch^{6,13}, Stany Chretien⁶, Luba Tcherpanov^{11,12}, Véronique Baudin-Creuzat¹⁴, Renaud Seigneuric⁸, Michaela Fontenay^{4,15}, Carmen Garrido^{7,8,16}, Olivier Hermine^{1,3,4,17} & Geneviève Courtois^{1,3,4}

PROTE

Cancer Research

VIROLOGY JOURNAL

JIM

PNAS

Cell Host & Microbe

nature immunology

PALL Life Sciences



Very strong Industry acceptance Pharma and Biotech





Molecular interactions in Biopharma



Discovery & Basic research

- Specificity = Which molecule?
 - Targets
 - Signalling mechanisms
 - Affinity = How strong is the binding?
 - Kinetics = How fast if will act?
 - Mapping = Where it will bind?

Bioprocess development / recombinant proteins

- Concentration = Do I have enough?
 - Clones, cell cultures conditions
 - Purity = Safe to use?
 - Formulation (Stability & Activity)


Clinical development

- Action or Immune reaction
 - ADA⁹




Real time approach : More Infomative

	Association	Dissociation	Ko = Affinty
Example A	1.00E+04	1.00E-02	1.00E-06
Example B	1.00E+02	1.00E-04	1.00E-06



Traditional routes to answers: with label and at equilibrium


Protein Quantitation with ELISA:






Total assay time: 3-4 hours plus overnight incubation

Requires 100-200 μ l sample

Target-specific antibodies



Capabilities of BLI Analytical Instruments




- Protein Concentration
 - Association rate
 - Dissociation rate
 - Binding Affinity

Myszka et. al. *Direct Comparison of Binding Equilibrium, Thermodynamic and Rate Constants Determined by Surface- and Solution-based Biophysical Methods*
Protein Science 2002, 11 (5), 1017-1025.




Octet Versatility in Interaction Analysis



13




What is Interferometry?




14



Measuring of interactions between waves





Two waves in phase
→ Constructive interference



Two waves 180° out of phase
→ Destructive Interference

Bio-Layer Interferometry (BLI)


- Much studied physical principle commonly used in surface profiling, semiconductor industry, astrophysics
- Optical layer reflects simple white light; second reflection from tip of biosensor, both reach detector
- Analyte binding changes thickness of bio-layer, which is measured at detector




16



Monitoring the Interference: Real Time






$$\text{Intensity } \lambda = f(\lambda, \ell)$$



17



Capabilities of BLI Analytical Instruments




- Protein Concentration
- Association rate
- Dissociation rate
- Binding Affinity

Myszka et. al. *Direct Comparison of Binding Equilibrium, Thermodynamic and Rate Constants Determined by Surface- and Solution-based Biophysical Methods*
Protein Science 2002, 11 (5), 1017-1025.


19



Label-free Techniques



- SPR (Surface Plasmon Resonance)
- QCM (Quartz Crystal Microbalance)
- SAW (Surface Acoustique Wave)





- They differ in detection method, but they all provide « Sensograms »
- They largely differ in speed, throughput, sample compatibility, usability

18



The BLI Sensor



Smaller ligand size allows greater packing on surface – greater binding capacity

20



Current Biosensors

Generic Biosensors:

- Streptavidin (SA)
- Super Streptavidin (SSA)
- High Precision Streptavidin (SAX)
- Amine Reactive 2nd Generation (AR2G)
- Aminopropylsilane (APS)

Tag Biosensors:

- Anti-Penta-HIS (HIS)
- Ni-NTA (HIS)
- Anti-GST
- Anti-FLAG

Current Biosensor Kits:

- Immunogenicity (ADA)
- Residual Protein A
- CHO HCP

Antibody Biosensors:

- Anti-Human IgG Fc (AHQ)
- Anti-hIgG Fc Capture Surface (AHC)
- Anti-Murine IgG Fv (AMQ)
- Anti-Murine IgG Fc Capture (AMC)
- Protein A
- Protein G
- Protein L
- Anti human IgG Fab




Custom Sensors

- Custom Capture Sensor
- Custom HCP kit


21



Sensor regeneration




BLI: One Platform, Many Capabilities



Quantitation

- Direct, 1-step
- Sandwich
- ELISA
- mg/mL – pg/mL

Kinetics

- k_a , k_d , K_D
- Proteins
- Peptides, Oligos
- Small molecules

Specificity

- Function testing
- Rank ordering
- Epitope binning
- Isotyping

23



BLI solutions A Wide range of applications

Kinetic Applications

- Protein - protein interactions
- Protein - small molecule interactions
- Liposome - protein/antibody binding
- Bacteria – antibody binding
- Virus-like particle - protein binding
- DNA aptamer binding
- Glycan – protein binding
- GPCR-Protein binding

Quantitative Applications

- Titer determination
- Rapid protein IgG quantitation
- Quantitation assays for ELISA replacement
- Residual Protein A contamination
- Protein/Antibody Quantitation
- Plant protein quantitation in crude extracts
- Host-cell protein contamination
- Immunogenicity (low and high affinity ADA's)

Screening Applications


- Screening proteins for crystallization studies
- DNA aptamer screening
- Small molecule fragment screening
- Secondary screening and hit validation
- DNA-DNA mismatch detection
- Phage binding (phage display)
- Protein/peptide/small molecule inhibition
- Clone selection in media
- Monitoring protein expression
- Bioreactor monitoring
- Epitope mapping/binning

Assay Development Applications

- Media development
- Process development
- Antibody subtyping
- Antibody pair selection



FortéBio OCTET application field : Innovator




25

PALL Life Sciences



FortéBio OCTET application field : Biosimilar



Historical applications

Second wave

Most Recent



Pall ForteBio is a Leader in Label-Free Protein Analysis



- Full life-cycle offering for protein and other biomolecular analysis
 - Instruments, consumables, software, service
 - Label-free assays based on Biolayer Interferometry (BLI)

27

PALL Life Sciences



Octet Systems



Analyte MW: > 150 D	Analyte MW: > 5 kD on Octet QK ^e , QK384	
LOD: 25 ng/mL for IgG on Protein A biosensors in 5 min assay	LOD: 100 ng/mL for IgG on Protein A biosensors in 5 min assay	
		<ul style="list-style-type: none"> ▪ 16-channel ▪ 96-, 384-well ▪ ≥ 40 µL sample volume ▪ 2 plate positions ▪ Biosensor re-racking ▪ Automation-friendly
		<ul style="list-style-type: none"> ▪ 8-channel ▪ 96-well ▪ 1 plate position ▪ Biosensor re-racking*
		*Not available on Octet QK ^e

28

PALL Life Sciences



The Classic Octet System




- Optics box moves sensors to samples
- Samples stay inside microplate; non-destructive testing
- “Dip and Read” format allows a large number of interactions to be studied in one experiment

29

Life Sciences



Unattended throughput



supports 96 sensors two 96-well or 384-well sample plates

30

Life Sciences



Octet Systems : K2 – Two channel automated




- 2 Channel Read Head
- 96 Well Microplates
- >180 µL Assay Volume
- 1 Microplate Position
- Biosensor Reracking

31

Life Sciences



What Is the Octet K2 System?



32

Life Sciences



Octet HTX – Unmatched Performance and Speed



Flagship of the Octet Line

Epitope Binning: 2-8 x Faster

Full Plate Quantitation in 2 Min


Full Plate Kinetic Screen in Mins, Not Hrs

33

PALL Life Sciences



Breakthrough Technology: User-Selectable Read Head



34



Automated Use of Octet 384 Systems



- Clone selection and screening
- Epitope mapping and binning
- Off-rate screening for antibody ranking
- Host cell protein and residual protein A detection
- Immunogenicity monitoring
- Small molecule and fragment screening

35

PALL Life Sciences



BLItz : The personal Analyzer



Label-free Assays in a Drop



36

PALL Life Sciences



BLITZ: Brilliantly Simple



3 Simple Steps:


- 1. Pipet a drop**
- 2. Insert a biosensor**
- 3. Close cover to start assay**

37

PALL Life Sciences



Workflow for Quantitation




- The binding rates of test samples are measured and interpolated from the standard curve to determine concentration
- 96 samples analyzed in 15 - 30 minutes
- Reuse of standard curve is optional

39

PALL Life Sciences



BLI: One Platform, Many Capabilities



Quantitation

- Direct, 1-step
- Sandwich
- ELISA
- mg/mL – pg/mL

Kinetics

- k_a , k_d , K_D
- Proteins
- Peptides, Oligos
- Small molecules

Specificity

- Function testing
- Rank ordering
- Epitope binning
- Isotyping



38

PALL Life Sciences



Octet Advantages in Bioprocessing Applications

Method Assay Steps



Octet unaffected by crude sample


IgG expressed in E. Coli

Octet Centrifuged	Octet Crude
0 $\mu\text{g/ml}$	0 $\mu\text{g/ml}$
280 $\mu\text{g/ml}$	281 $\mu\text{g/ml}$
42 $\mu\text{g/ml}$	47 $\mu\text{g/ml}$
323 $\mu\text{g/ml}$	327 $\mu\text{g/ml}$
132 $\mu\text{g/ml}$	130 $\mu\text{g/ml}$

40

PALL Life Sciences

Direct quantification : using catalogue sensors




Custom quantification


- A custom quantification consisted in creating a sensor (from one in the sensor catalog) with a surface able to bind the molecule of interest in the media
 - This quantification can be split into two phases (as shown in the schema)
 - Loading of a ligand (often antibody) able to bind the analyte we aim to measure
 - Immunocapture (HIS TAG, AMC, AHC, ...)
 - “Covalent” loading (SA, AR2G)
 - Dipping into the analyte solution and measure of the association

Ex: basic protocol

Use catalog sensor




Detection of 24kDa Viral Protein on Octet QKe



Immobilized capture antibody onto Streptavidin Biosensors and looked at binding of sample for ten minutes

Able to detect binding of protein down to 24 ng/mL with one step quantitation assay


Quantitation of HIS-tagged PAI (43 kDa) on Octet RED



Rapidly quantitate His-tagged proteins on the Octet using Anti-Penta HIS Biosensors



Octet Speed & Reproducibility Comparison to ELISA



Octet = 30 minutes,
ELISA = 5-6 hours

Typical RSD values for
multiple users and
systems:


ELISA 5 - 20%
Octet < 5%

Data presented at IBC Antibody Production, 2008 by Keith Davis from Pfizer MO.

45



Expression Ranking : what are the best conditions



Rank	Binding Rate	%CV
Sample 1	1.991	10.8%
Sample 7	1.646	2.4%
Sample 3	1.373	3.5%
Sample 9	0.864	3.2%
Sample 11	0.614	2.2%
Sample 6	0.230	3.7%
Sample 5	0.112	0.6%
Sample 12	0.031	4.3%
Sample 8	0.015	2.0%
Sample 10	0.004	4.3%
Sample 4	0.002	8.9%
Sample 2	0.000	NA

Highest expression
↑
Lowest expression

47



Case Example : Bio Process development



+




Develop the cultures

Analyze the results

46



BLI: One Platform, Many Capabilities



Quantitation

- Direct, 1-step
- Sandwich
- ELISA
- mg/mL – pg/mL

Kinetics

- k_a , k_d , K_D
- Proteins
- Peptides, Oligos
- Small molecules

Specificity


- Function testing
- Rank ordering
- Epitope binning
- Isotyping

48






An Ideal Sensogram




In a simple 1:1 binding model, the association and dissociation phases are described by a single exponential function

49




Putting it together : from sensogram to K_D

For a simple 1:1 Interaction:




Using successive iterations until a fit is found :




k_a = rate of association or "on-rate"

k_d = rate of disassociation or "off-rate"


51




An Ideal Sensogram



50




Kinetics Work Flow



52



Binding Kinetics: Effortless



Application areas:

- Protein-protein interactions
- Protein quality through binding assays
- Biotherapeutic drug development

Analysis methods complemented/replaced:


- SPR : faster, simpler, more flexible
- ELISA : faster, more accurate
- ITC: faster, different insight

53



Parathyroid Hormone (PTH) Binding to PTH Receptor

(Octet) Biotin-PTH1R ECD on SA biosensor – PTH(1-34)NH₂ in solution.
(Biacore) PTH1R on CM5 chip – PTH(1-34)NH₂ in solution.



Instrument	K_D
Biacore X	$4.9 \mu\text{M}$
ForteBio Octet RED	$2.8 \mu\text{M}$
Microcal iTC	$3.4 \mu\text{M}$

PTH is of clinical interest for its ability to stimulate bone formation.

(Octet) *J. Biol. Chem.* **2009**, *284*, 28382-28391.


(Biacore) *Biochemistry*, **2000**, *39*, 8878-8887.



DNA-Binding Protein Kinetic Analysis on the Octet QK




DNA-Binding Protein Kinetic Analysis on the Octet QK



Binding of DNA-Binding Protein to Immobilized Biotinylated ssDNA

Rapid Analysis of Binding of Transcription Factors and Other Promoter Elements to Specific DNA Sequences

55



Binding of DNA-Binding Protein to Immobilized Biotinylated ssDNA


Rapid Analysis of Binding of Transcription Factors and Other Promoter Elements to Specific DNA Sequences

56





Small molecule screening/binding characterization



Data very consistent across both platforms


* Small Fragments binding to Kinase- courtesy of Roche

57

PALL Life Sciences



Determination of the KD for Ab – Ag interactions (K2 in SWU)



59

PALL Life Sciences



Round Robin on SPR & BLI

David Myszka Publication Biosensor Benchmark Study Comparison of Biacore Models vs. Octet QK

50 kDa Fab binding to a 60 kDa antigen				
Instrument	Model	n	k_a (M-1 s-1)	k_d (M-1 s-1)
Biacore	A100	3	(1.1 ± 0.3) X 10^5	(0.60 ± 0.26) X 10^4
	T100	33	(1.4 ± 0.8) X 10^5	(0.43 ± 0.14) X 10^4
	S51	8	(1.7 ± 1.4) X 10^5	(0.86 ± 0.75) X 10^4
	3000	77	(1.3 ± 0.9) X 10^5	(0.75 ± 1.2) X 10^4
	2000	76	(1.4 ± 1.7) X 10^5	(0.57 ± 1.0) X 10^4
	1000	2	(0.82 ± 0.19) X 10^5	(1.3 ± 0.7) X 10^4
	X100	2	(1.26 ± 0.08) X 10^5	(0.38 ± 0.03) X 10^4
	X	8	(0.87 ± 0.29) X 10^5	(0.65 ± 0.20) X 10^4
	Flexchip	4	(1.2 ± 0.2) X 10^5	(0.46 ± 0.10) X 10^4
ForteBio	Octet QK	3	(1.7 ± 0.6) X 10^5	(0.70 ± 0.22) X 10^4
Study Average		258	(1.4 ± 1.3) X 10^5	(0.61 ± 0.87) X 10^4
				0.62 ± 0.98

"These results demonstrate that when this biosensor assay was designed and executed appropriately, the reported rate constants were consistent, and independent of which protein was immobilized and which biosensor was used."

Riche et. al. Analytical Biochemistry 386 (2009) 194–216

Fast. Accurate. EASY.

58


fortéBIO
A Division of Pall Life Sciences



Determination of the KD for Ab – Ag interactions

Project 2

▪ Loading at 5ug/ml mAb



Good fitting (global 1:1). By excluding the highest concentrations, it gets even slightly better.


Sample ID	Conc. (nM)	KD (M)	KD Error	k_{on} (1/Ms)	k_{on} Error	k_{dis} (1/s)	k_{dis} Error	Rmax	Rmax Error	k_{obs} (1/s)	Req	Req/Rmax(%)	Full X^2	Full R^2
	200	5.25E-10	<1.0E-12	1.34E+05	1.61E+02	7.02E-05	4.82E-08	0.4338	0.0001	2.68E-02	0.4327	99.7	0.402971	0.99727
	100	5.25E-10	<1.0E-12	1.34E+05	1.61E+02	7.02E-05	4.82E-08	0.4308	0.0001	1.34E-02	0.4285	99.5	0.402971	0.99727
PG VB26	50	5.25E-10	<1.0E-12	1.34E+05	1.61E+02	7.02E-05	4.82E-08	0.41	0.0001	6.76E-03	0.4058	99	0.402971	0.99727
	25	5.25E-10	<1.0E-12	1.34E+05	1.61E+02	7.02E-05	4.82E-08	0.4199	0.0002	3.41E-03	0.4113	98	0.402971	0.99727
	12.5	5.25E-10	<1.0E-12	1.34E+05	1.61E+02	7.02E-05	4.82E-08	0.4199	0.0003	1.74E-03	0.403	96	0.402971	0.99727

PALL Life Sciences



Determination of the KD for Ab – Ag interactions

- Loading at 2.5ug/ml mAb



Good fitting (global 1:1). By excluding the highest concentrations, it gets even slightly better.


Sample ID	Conc. (nM)	KD (M)	KD Error	kon(1/Ms)	kon Error	kdis(1/s)	kdis Error	Rmax	Rmax Error	kobs(1/s)	Req	t _{eq} /Rmax(%)	Full X ²	Full R ²
PG VB26	200	5.63E-10	1.23E-12	1.30E+05	2.48E+02	7.31E-05	7.80E-08	0.3077	0.0001	2.60E-02	0.3069	99.7	0.489537	0.995908
	100	5.63E-10	1.23E-12	1.30E+05	2.48E+02	7.31E-05	7.80E-08	0.3123	0.0001	1.31E-02	0.3105	99.4	0.489537	0.995908
	50	5.63E-10	1.23E-12	1.30E+05	2.48E+02	7.31E-05	7.80E-08	0.2881	0.0001	6.57E-03	0.2849	98.9	0.489537	0.995908
	25	5.63E-10	1.23E-12	1.30E+05	2.48E+02	7.31E-05	7.80E-08	0.2701	0.0002	3.32E-03	0.2641	97.8	0.489537	0.995908
	12.5	5.63E-10	1.23E-12	1.30E+05	2.48E+02	7.31E-05	7.80E-08	0.2177	0.0003	1.70E-03	0.2083	95.7	0.489537	0.995908

61



Determination of the KD for Ab – Ag interactions

- Loading at 1ug/ml mAb



Good fitting (global 1:1). By excluding the highest concentrations, it gets even slightly better.

Sample ID	Conc. (nM)	KD (M)	KD Error	kon(1/Ms)	kon Error	kdis(1/s)	kdis Error	Rmax	Rmax Error	kobs(1/s)	Req	Req/Rmax(%)	Full X ²	Full R ²
PG VB26	200	7.86E-10	2.64E-12	1.50E+05	4.72E+02	1.18E-04	1.33E-07	0.1568	0.0001	3.00E-02	0.1562	99.6	0.242063	0.994755
	100	7.86E-10	2.64E-12	1.50E+05	4.72E+02	1.18E-04	1.33E-07	0.1644	0.0001	1.51E-02	0.1632	99.3	0.242063	0.994755
	50	7.86E-10	2.64E-12	1.50E+05	4.72E+02	1.18E-04	1.33E-07	0.0765	0.0001	7.60E-03	0.0753	98.4	0.242063	0.994755
	25	7.86E-10	2.64E-12	1.50E+05	4.72E+02	1.18E-04	1.33E-07	0.1228	0.0001	3.86E-03	0.119	96.9	0.242063	0.994755
	12.5	7.86E-10	2.64E-12	1.50E+05	4.72E+02	1.18E-04	1.33E-07	0.0648	0.0001	1.99E-03	0.061	94.1	0.242063	0.994755


62



Affinity Characterization of IgG Products (8 channel Octet)

96 supes can be assayed in 10-15 minutes with the HTX

16 mutant culture supernatants were assayed in triplicate on Ag-Fc fusion protein bound to Protein A biosensors in less than 30 minutes



These crude samples were unsuitable for Biacore analysis due to media composition and unfilterable cell debris.


Rapidly eliminated 5 out of 16 mutants by off rate criteria.

Data from Brian Miller at Biogen IDEC

63



Binding of mAbs in PBS to VLPs on Octet RED



* Data courtesy of Integral Molecular

64






Kinetics of *Arabidopsis* Proteins measured on Octet QK

Autophosphorylation of Tyr-610 in the receptor kinase BAK1 plays a role in brassinosteroid signaling and basal defense gene expression

Man-Ho Oh^{a,b}, Xiaofeng Wang^a, Xia Wu^a, Youfu Zhao^a, Steven D. Clouse^a, and Steven C. Huber^{a,b,1}

^aAgricultural Research Service, United States Department of Agriculture and ^bDepartment of Plant Biology, University of Illinois, Urbana, IL 61801; ¹Department of Horticultural Science, North Carolina State University, Raleigh, NC 27695; ²Physiological and Molecular Plant Biology Program, Department of Plant Biology, University of Illinois, Urbana, IL 61801; and ³Department of Crop Sciences, University of Illinois, Urbana, IL 61801

Edited by Jianming Li, University of Michigan, Ann Arbor, MI and accepted by the Editorial Board September 1, 2010 (received for review January 5, 2010)



Despite gross changes in plant growth rate, Octet QK data demonstrates the kinetic parameters of the mutated kinase are unchanged from wild type.

Octet QK parameters:


- Standard 96W plate
- Anti-Murine biosensors
- 30C, 1000 RPM in MOPS buffer
- Loaded anti-GST antibody to capture GST-BAK protein
- 15 minute association w/ BRI1
- 15 minute dissociation

PALL Life Sciences




Galectin 8 vs. VHH-H4 : Assay setup & Raw data

The following assay parameters were used for all other VHH studied vs. Galectin 8




Galectin 8 vs. VHH-H4 : Reference substracted data



PALL Life Sciences




Galectin 8 vs. VHH-H4 : KD calculation



Sample ID	Conc. (nM)	KD (M)	kon(1/Ms)	kdis(1/s)
VHH-H4	500	2.27E-07	5.42E+04	1.23E-02
VHH-H4	250	2.27E-07	5.42E+04	1.23E-02
VHH-H4	83	2.27E-07	5.42E+04	1.23E-02



BLI: One Platform, Many Capabilities



Quantitation

- Direct, 1-step
- Sandwich
- ELISA
- mg/mL – pg/mL

Kinetics

- k_a , k_d , K_D
- Proteins
- Peptides, Oligos
- Small molecules

Specificity

- Function testing
- Rank ordering
- Epitope binning
- Isotyping

69

Life Sciences



Assay Development: Smarter

Tests performed:

- Binding pair optimization
- Immunoassay development
- Cross-blocking studies

Analysis methods complemented/replaced:


- ELISA
- SPR

Application areas:

- Immunoassay development
- Academic research ELISAs
- Biomarker discovery and development
- Reagent quality analysis
- Antibody manufacturing
- Enzyme and other protein manufacturing


70

Life Sciences




Use of the Octet platform in protein crystallography

Data courtesy of Jiamin Yu from Novartis | CHI The Bioprocess Summit 2011 | ForteBio Workshop



- Aggregated and mis-folded proteins bind slower to a biosensor surface due to larger size

- Select for monodisperse proteins for crystallography based on a simple label-free binding assay using the Octet platform to do a Binding Rate Screen (BRS)



Monodisperse
Faster binding rate

Aggregated or Misfolded
Slower binding rate

Binding Rate Screen - a high-throughput assay in soluble lysate for prioritizing protein expression constructs
Anal Biochem. 2010 Apr 15;399(2):276-83.




Binding Rate Screen Identifies Constructs That Produce Crystals



Data courtesy of Jiamin Yu from Novartis | CHI The Bioprocess Summit 2011 | ForteBio Workshop

- SEC on crude lysate
- Analyze fractions with anti-his immuno assay
- Compare to Octet binding rate assay
- Constructs E & F crystallized

1 minute Octet crude lysate assay correctly prioritizes protein constructs for crystallization



Binding Rate Assay



7


71

Life Sciences


Life Sciences



Antibody Epitope binning/pairing



Subtyping human IgG's



Recover for MS : EPO quant trouble shooting

Customer :

Research at CFB

The Novo Nordisk Foundation Center for Biosustainability (CFB) at the Technical University of Denmark aims at developing new knowledge and technologies to help facilitate the transformation from the existing oil-based chemical industry to a more sustainable bio-based society, in which chemicals are produced biologically. Furthermore, CFB aims at accelerating the development of genome-scale science for CHO cell lines for improved bioprocessing.

The center is organized in 11 research sections and a Core which aims at platform strain and technology development purposes:

- Yeast Cell Factories (YCF)
- CHO Cell Line Engineering and Design (CLED)
- Glyco-Engineering of CHO (GEC)
- Network Reconstructions and in silico Biology (NRISB)
- Genome-Scale CHO in silico models (GSCiM)
- Yeast Synthetic Biology (YSB)
- High-throughput Molecular Bioscience (HTMB)
- Bacterial Cell Factories (BCF)
- New Bioactive Compounds (NBC)
- iLoop
- Cho-Core


Application : Customer has developed a custom quant assay for EPO (supernatants / purified)

Issue: in certain eukaryote supernatants the sensitivity is unexpectedly low compared to others, with no apparent reason. Teams suspects some issues with (non)specific binding of competing proteins

Solution : Recovery of capture and MS analysis




Recovery for MS : ligand fising






LC - MS

1. Neutralize recovery solution
2. Micro UPLC
3. MS-MS analysis



77



PALL Life Sciences



Octet Platform

fortéBIO

A Division of **Pall Life Sciences**

Top 10 Applications

- 1.Kinetic Characterization
- 2.Ab Titer
- 3.Cross-blocking and Ab Pairing (Epitopes)
- 4.Non-Ab Titer
- 5.Small Molecule**
- 6.Mechanism of Action
- 7.Reagent Development (Immunoassay)
- 8.Lipids, Membranes
- 9.Vaccines
- 10.Oligonucleotides


Source: North American customer survey of 478 application responses

78

PALL Life Sciences



High-Throughput Small Molecule Fragment Screening and Hit Validation on Octet RED384



- 16-channels
- 384 tilted-well plate (40µL min volume)
- 2 plate positions



79

PALL Life Sciences



Full Robotics Compatibility for True High-Throughput Label-Free Measurements




Integration-Ready for All Major Robot Vendors

80

PALL Life Sciences



Large Fragment Libraries Screened & Analyzed Quickly and Easily




81

PALL Life Sciences



Small Molecule Fragment Screening on Octet RED384 vs Biacore Small Fragments binding to Kinase




Data very consistent across both platforms

83

PALL Life Sciences



Hit Validation by Titration Series




82

PALL Life Sciences




Small Molecule Kinetics on Octet RED384 vs Biacore Roche Validation Study

Octet RED384



Biacore



238 Dalton fragment binding to biotinylated kinase;
6 point dilution series in TBS, 1mM DTT, 0.005%P20, 5%DMSO

Data Courtesy of Frank Podlaski, Discovery Technologies, Roche, Inc.

84

PALL Life Sciences




Confirmation of hits by X-Ray

Crystal structures – Orthogonal Method



- Co-crystallisation of top hits with target protein



Benefits of the Octet Platform

▪ Higher throughput, faster workflow

- Concentration measurements faster than ELISA or HPLC
- Screen more samples in less time than SPR
- Little dilution or **sample preparation required**
- **Full sample recovery**
- Fixed volume allows long on/off rate determinations

▪ Unparalleled ease-of-use

- Multiple users are easily trained to operate the system
- **Microfluidics-free** format and disposable biosensors
- **No maintenance**, cleaning, or instrument prep time required

▪ Versatility

- **Crude samples** can easily be measured
- **Wider solvent tolerance** and assay flexibility

▪ Price

- **Lower capital acquisition cost** than traditional SPR systems
- Lower usage/maintenance costs SPR systems
- **Lower cost per sample** than Protein A HPLC or commercial ELISA kits

86

PALL Life Sciences



Life Sciences

Better Lives. Better Planet.SM

ENABLING A
GREENER
FUTURESM

www.pall.com/green