

X-ray crystallography practical

Oleron 2017

Last update: 19 June 2017

Table of Contents

Schedule of the X-ray crystallography practical.....	1
Data collection on a Gd derivative on beamline FIP.....	2
Data reduction with iMosflm.....	3
Start iMosflm.....	3
Indexation.....	3
Integration and scaling.....	5
Data reduction with XDS package.....	6
Introduction.....	6
Default xds parameter files.....	8
Run all step in a row.....	8
Re-run scaling.....	8
xscale and xdsconv.....	9
Data reduction with XDS in graphic/automated mode.....	11
Graphic mode, using xdsGUI.....	11
Automated data reduction.....	12
Quick SAD phasing with Phaser in ccp4i.....	14
Quick SAD phasing in ccp4i2.....	18
Run the SHELX automated pipeline.....	18
Results analysis.....	20
Automated SAD phasing with Phaser.....	23
SAD phasing, step by step, with ccp4i.....	25
MR.....	34
Complete model.....	34
Partial/remote model.....	34
Refinement / construction.....	34
Complete model.....	34
Partial/remote model.....	34

Schedule of the X-ray crystallography practical

Tuesday 20 June, 9:15AM-12:45AM (1h45 + 1h30)

9h00-9h15 : sponsors

Présentation du TP, constitution des groupes (~15 min)

Data reduction (~30 min)

Données : BM30A-2014-11-19 (native, 164 frames, 1.1 Å) lyso poule
15-05-13-lyso/lyso-Gd_SAD (anomales Gd, 300 frames, 1.65 Å)
Soft : xds, xdsGUI ou mosflm
Support : tutorial_xds ou tutorial_iMOSFLM

Phasage par MR (~30 min)

Données : BM30A-2014-11-19 , divers modèles : lyso incomplet de poule ou
modèle complet "basse" homologie (boeuf)
Soft : phaser
Support : ?

Phasage SAD (~30 min)

Données : lyso-Gd_SAD
Soft : ccp4
Support : tutorial_SAD, tutorial_SAD_bis

Pause café (11h00-11h15)

Construction des parties manquantes / affinement (~45 min)

Données : BM30A-2014-11-19, modèle lyso incomplet
Soft : coot, refmac
Support : ?


Construction des parties manquantes / affinement (~45 min)

Données : lyso-Gd_SAD, modèle lyso incomplet, modèle ligand+Gd
Soft : coot, refmac
Support : tutorial_SAD, tutorial_SAD_bis


Repas (12h45)

Data collection on a Gd derivative on beamline FIP

The fluorescence of Gd was measured with a Roentec MCA at the Gd LIII edge. Raw data are in
Edge/lyso_1_Gd1 (columns 5 and 7)
and the plot vs Energy in
Edge/lyso_1_Gd1.jpg



The spectrum was processed with Chooch. Final drawing of calculated f' and f'' is in
Edge/final.jpg



Based on that, beam energy was tuned to 7242.6 eV, and a single-wavelength dataset was collected (300 frames, 1 deg each). Frames (compressed with bzip2) are named
img/e000_prefix_1_00xxx.img.bz2

Select the correct space-group: nb 96, $P4_32_12$

Session Settings Help

157.59 157.80 114.20 5.00 10.0 0.92 0.71 0.00 20 292

Autoindexing

e000_prefix_1_####.img: 1-2, 89-91

Image	ε range	Auto	Man	Del	> I/ε (I)	Find	Use
1	0.00 - 1.00	773	0	0	673	☒	☒
2	1.00 - 2.00	880	0	0	745	☒	☒
89	88.00 - 89.00	818	0	0	687	☒	☒
90	89.00 - 90.00	784	0	0	658	☒	☒
91	90.00 - 91.00	733	0	0	643	☒	☒
Total		3988	0	0	3406		

Lattice 1

Solution	Lat.	Pen.	a	b	c	ε±	ε'	εz	ε(x,y)	Nref	ex beam
1 (ref)	aP	0	37.8	75.6	75.7	89.7	89.8	89.7	0.73	3268	0.47 (0.2)
2 (ref)	aP	1	37.8	75.6	75.7	90.3	89.8	90.3	0.73	3268	0.47 (0.2)
3 (ref)	mP	2	37.8	75.8	75.5	90.0	90.4	90.0	0.71	3238	0.44 (0.2)
4 (ref)	mC	2	107.3	106.8	37.7	90.0	90.5	90.0	0.86	3306	0.38 (0.2)
5 (ref)	mP	3	37.8	75.6	75.7	90.0	89.7	90.0	0.75	3257	0.48 (0.2)
6 (ref)	mP	3	75.7	37.8	75.8	90.0	90.5	90.0	0.83	3305	0.52 (0.2)
7 (ref)	mC	4	107.3	106.8	37.7	90.0	90.5	90.0	0.86	3306	0.38 (0.2)
8 (ref)	oP	4	37.8	75.6	75.8	90.0	90.0	90.0	0.76	3249	0.49 (0.2)
9 (ref)	oC	4	106.8	107.7	37.8	90.0	90.0	90.0	0.99	3335	0.33 (0.2)
10 (ref)	tP	5	75.7	75.7	37.8	90.0	90.0	90.0	0.77	3263	0.48 (0.2)
11 (reg)	mC	110	155.5	37.7	75.7	90.0	90.1	90.0	-	-	-
12 (reg)	mC	110	37.7	155.5	75.7	90.0	90.2	90.0	-	-	-
13 (reg)	oC	111	37.7	155.5	75.7	90.0	90.0	90.0	-	-	-
14 (reg)	mC	111	37.7	155.8	75.5	90.0	90.4	90.0	-	-	-

Lattices:

Spacegroup: P43212 Prior cell:

Mosaicity:

Green warnings 0

then run “Cell Refinement”

Session Settings Help

157.59 157.80 114.20 5.00 10.0 0.92 0.71 0.00 20 292

Cell refinement

Lattice: e000_prefix_1_####.img: 1-4, 91-94

Parameter	Value	Fix
Beam x	157.55	☐
Beam y	157.87	☐
Distance	114.20	☐
Y-scale	1.0003	☐
Tilt	-0.22	☐
Twist	0.15	☐
Tangential offset	0.000	☐
Radial offset	0.000	☐
RMS residual	0.092	☐
RMS res. (central)	0.066	☐
RMS res. (weighted)	0.040	☐

RMS residual

Parameter	Value	Fix
ε(x)	-0.02	☐
ε(y)	-0.05	☐
ε(z)	-0.04	☐
a	77.27	☐
b	77.27	☐
c	38.70	☐
ε±	90.00	☐
ε'	90.00	☐
εz	90.00	☐
Mosaicity	0.404	☐


Cell

Cell	a	b	c	ε±	ε'	εz
Initial	75.70	75.70	37.79	90.00	90.00	90.00
Final	77.27	77.27	38.70	90.00	90.00	90.00
Std dev	0.00	0.00	0.00	0.00	0.00	0.00

Green warnings 0

Integration and scaling

Select “Integration”



Then run “QuickSymm” (quick scaling to check the symmetry) and “QuickScale” (scaling)
Check the Rmerge and Rmeas

Notably, compare Rmerge or Rmeas calculated either on:
all I+ and I- (includes the centrosymmetry (Fridel's law))
or
within I+ or I- (excludes the centrosymmetry (Fridel's law))

DeltaAnom correlation between half-sets tells you whether the anomalous signal is significant: You divide your data into two sets, and compare the anomalous difference I+ - I- in the two sets. If the correlation is high, the signal is significant.

You can look also at the other figures on the table. All the integrated and scaled data are stored in the file **ctruncate_A5_1_00002-unique.mtz**
This file contains h, k, l, FreeRflag, F, SigF, Dano, SigDano, F(+), SigF(+), F(-), SigF(-), Isym, Imean, SigImean I(+), SigI(+), I(-), SigI(-)

For the following steps, you can use the **ctruncate_A5_1_00002-unique.mtz** file.

Data reduction with XDS package

Introduction

XDS is a suite of programs dedicated to the reduction of macromolecular crystallography data. The suite of programs include:

xds: data processing, from images to unmerged h,k,l Intensities, $\sigma(I)$

xscale: scaling and merging Intensities from either one or several data sets.

xdsconv: converts reflection data files as obtained from xds or xscale into various formats required by software packages for crystal structure determination like CCP4, CNS (X-PLOR), or SHELX.

2cbf: converts a detector image file to CBF format. (not often used)

merge2cbf: converts a series of detector image files to CBF format. (not often used)

cellparm: used to determine the mean of the cell parameters obtained from processing several data sets from the same crystal form. (not often used).

Only **xds**, **xscale** and **xdsconv** will be used here.

xds requires

- diffraction images
- a parameter file called XDS.INP that contains all the necessary information regarding the experimental setup.

Most of the time, an XDS.INP file is generated automatically when you launch a data collection at a synchrotron. However, the file XDS.INP needs some editing during the data processing, but only a few input parameters require to be looked at. See the commented XDS.INP file for further details.

The whole data processing includes 7 steps defined in the JOB= command line. Each step generates a log file named with the .LP suffix.

XYCORR: computes a table of spatial correction values for each pixel: allow to precisely localise each pixel of the detector. Fully automatic, to be done once.

files created:

X-CORRECTIONS.cbf

Y-CORRECTIONS.cbf

XYCORR.LP

INIT: determines an initial background for each detector pixel and finds the trusted region of the detector surface. Needs 5 to 10 images to run properly (look at BACKGROUND_RANGE= command). To be done once.

files created:

BKGINIT.cbf

BLANK.cbf

GAIN.cbf

INIT.LP

COLSPOT: collects strong diffraction spots from a specified subset of the data images (see

SPOT_RANGE= command).

files created:

FRAME.cbf
SPOT.XDS
COLSPOT.LP

IDXREF: interprets observed spots by a reciprocal lattice and refines all diffraction parameters (cell dimensions, orientation matrix, crystal-detector distance, etc ...).

files created:

XPARM.XDS
IDXREF.LP

DEFPIX: defines the trusted region of the detector, recognizes and removes shaded areas, and eliminates regions outside the resolution range defined by the user.

files created:

BKGPIX.cbf
ABS.cbf
DEFPIX.LP

XPLAN: helps planning data collection. Tells you what data to collect in order to get the most complete data set. Only useful when at the synchrotron beamline, before launching the data collection.

files created:

XPLAN.LP

INTEGRATE: collects 3-dimensional profiles of all reflections occurring in the data images and estimates their intensities

files created:

INTEGRATE.HKL
INTEGRATE.LP

CORRECT: corrects intensities for decay, absorption and variations of detector surface sensitivity, reports statistics of the collected data set and refines the diffraction parameters using all observed spots.

files created:

ABSORP.cbf
DECAY.cbf
DX-CORRECTIONS.cbf
DY-CORRECTIONS.cbf
GX-CORRECTIONS.cbf
GY-CORRECTIONS.cbf
MODPIX.cbf
GXPARM.XDS
XDS_ASCII.HKL
CORRECT.LP

The different steps are presented in a series of directories, for sake of clarity. In practice, they can be performed in a single directory by successive modification of the input files and running the XDS package programs at the command line.

Default xds parameter files

Go to the directory **xds_step0_default-XDSINP**

You will find there the initial parameter file for xds, as created automatically by the beamline control software (simplified version)

```
xds_step0/XDS.INP
```

as well as a fully commented one

```
xds_step0/XDS.INP_sav
```

Run all step in a row

Go to directory **xds_step1_all** (783 sec elapsed time on a HP ElitBook 840 i5 notebook, with bziped frames; 252 sec elapsed time with uncompressed frames)

As we expect anomalous signal, the Friedel mates will differ. So uncomment the line

```
FRIEDEL'S_LAW=FALSE
```

in XDS.INP. Then, run

```
xds
```

at the command line (or `xds_par` for the paralleled version).

Re-run scaling

Go to the directory **xds_step2_CORRECT** (60 sec elapsed time)

At the previous step, xds automatically figures out the Bravais symmetry, and picked up space group P422 as a representative. To check for extinctions (helices), just select space group 96 (P4(3)2(1)2) and enter refined cell parameters in XDS.INP

```
SPACE_GROUP_NUMBER=96
```

```
UNIT_CELL_CONSTANTS= 77.268 77.268 38.704 90.000 90.000 90.000
```

and select only the final scaling step of the processing (CORRECT)

```
JOB= CORRECT
```

Then run

```
xds
```

In CORRECT.LP, check for low intensity of reflections that should be absent (marked with "*") in the list above lines

```
AVERAGE INTENSITY FOR 207 REFLECTIONS WHICH SHOULD  
BE SYSTEMATICALLY ABSENT IS 0.2% OF MEAN INTENSITY
```

xscale and xdsconv

Go to the directory **xds_step3_XSCALE-XDSCONV** (5.6 sec elapsed time)

Optional: use XSCALE for final scaling, merging of several dataset. Create the XSCALE.INP input file with the following lines. Not necessary for a single dataset, apart to redefine boundaries of resolution shells.

```
OUTPUT_FILE=XSCALE.HKL
INPUT_FILE= XDS_ASCII.HKL
```

and run

```
xscale
```

at the command line (alternatively: xscale_par).

Use XDSCONV to generate reflection files in CCP4 FP/DANO format (F, SigF, Dano, SigDano):

```
INPUT_FILE=XSCALE.HKL XDS_ASCII
OUTPUT_FILE=temp_ccp4.hkl CCP4
FRIEDEL'S_LAW=FALSE
```

and run

```
xdsconv
```

at the command line.

XDSCONV generates the input file F2MTZ.INP needed by f2mtz (CCP4 package) for the final conversion to binary mtz format. To run the CCP4 programs f2mtz just type the command:

```
f2mtz HKLOUT temp_ccp4.mtz < F2MTZ.INP
```

Use XDSCONV again to generate reflection files in CCP4 F+/F- format (F, SigF, F+, SigF+, F-, SigF-):

```
INPUT_FILE=XSCALE.HKL XDS_ASCII
OUTPUT_FILE=temp_ccp4_f.hkl CCP4_F
FRIEDEL'S_LAW= FALSE
GENERATE_FRACTION_OF_TEST_REFLECTIONS=0.05
```

and run again

```
xdsconv
```

at the command line

Then run CCP4 program f2mtz

```
f2mtz HKLOUT temp_ccp4_f.mtz < F2MTZ.INP
```

Then, to run cad (to convert indices to the CCP4-asymmetric unit),

```
cad HKLIN1 temp_ccp4.mtz HKLIN2 temp_ccp4_f.mtz HKLOUT Lyso-
Gd_SAD.mtz <<EOF
```

```
LABIN FILE 1 E1=FP    E2=SIGFP    E3=DANO E4=SIGDANO E5=ISYM
LABIN FILE 2 E1=F(+)  E2=SIGF(+)  E3=F(-) E4=SIGF(-) E5=FreeRflag
END
EOF
```

Data reduction with XDS in graphic/automated mode


Graphic mode, using xdsGUI

In xds_xdsGUI


Run

xdsGUI

at the command line




From there, all steps described above can be performed, starting with the edition of the XDS.INP parameter file




and then running xds (the “Run XDS” button).

Log files are available in the “XYCORR”, “INIT”, ... windows, with graphical display of the

statistics:



extra statistics are available running xdsstat in the “statistics” window:



Access to xscale and xdsconv is provided through the corresponding windows.

Automated data reduction

A series of automated scripted layers are available for automated data reduction with xds, such as

xdsme

xdsapp

Just run from the top directory

```
xdsapp -a -cmd --fried=false
```

The dataset will be automatically detected and processed (parameters from the header of images). Alternatively, xdsapp is also a graphic interface for xds.

Xia

Quick SAD phasing with Phaser in ccp4i

Launch

`ccp4i`

and define a new project with

`lyso-Gd_SAD/ccp4_SAD`

as working directory

Select “Phaser SAD Pipeline” (button highlighted in blue in Figure 1) (~330 sec elapsed time).

Then enter the following parameters:

- reflexion file: `xds_step3_XSCALE-XDS CONV/Lyso-Gd_SAD.mtz` (field highlighted in blue in Figure 2)
- sequence in fasta format directory `ccp4_MR` (field highlighted in green in Figure 2)
- heavy atom type: `GD` (field highlighted in red in Figure 2)
- wavelength: `1.7119` (field highlighted in orange in Figure 2)




Figure 1

Maximum Likelihood Experimental Phasing Initial parameters from /home/data/15-05-13-lyso/lyso-Gd_S...

Job title [No title given]

Mode for experimental phasing Single-wavelength anomalous dispersion (SAD)

Phaser SAD pipeline

Run SHELXD before Phaser

Run Parrot (density modification) after Phaser

Run Buccaneer (model building) after Parrot

Define data

MTZ in Full path.. /home/data/15-05-13-lyso/lyso-Gd_SAD/xds_step3/Lyso-Gd

Crystal unknown belonging to Project unknown

Dataset name unknown140515

F(+) F(+) SIGF(+) SIGF(+)

F(-) F(-) SIGF(-) SIGF(-)

FP FP SIGFP SIGFP

FREER FreeRflag

Resolution 38.635 A to 1.648 A; Wavelength 1.7119

Space group read from mtz file 'P 43 21 2'; Enantiomorph choice Both

Enter scattering from fluorescence scan (default is to calculate f' and f'' from wavelength)

Define atoms

Atom sites run ShelxC/D

Find 4 heavy atoms of type GD

Crystal contains cluster compound

LLG-map completion on: all atom types

Shelx parameters

Composition of the asymmetric unit

Total scattering determined by components in asymmetric unit

Component #1 protein sequence file Number in asymmetric unit 1

SEQ file Full path.. /home/data/15-05-13-lyso/lyso-Gd_SAD/ccp4_MR/193L.fast

Edit list Define another component

Accessory parameters

Output control

Expert parameters

Run Save or Restore Close

Figure 2

=> ~80% of residues built automatically

Upon completion of the job, and to analyze the log file, select the “Phaser_EP” job in the list (button highlighted in blue in Figure 3). Then, from the “View Files from Job”, select “View Job Results (new style)” (button highlighted in red in Figure 3)




Figure 3

Experimental map, sub-structure of anomalous atoms and model can be displayed with Coot:

Run

coot

at the command line, and load pdb files (button highlighted in blue in Figure 4) and mtz files (button highlighted in red in Figure 4) as listed below:




Figure 4

Sub structure of Gd atoms is in

lyso-Gd_SAD/LysoGd_SAD_3.1.pdb

Experimental map is in

lyso-Gd_SAD/LysoGd_SAD_3.1.mtz

and after automated density modification with parrot

lyso-Gd_SAD/LysoGd_SAD_3_parrot.mtz

The model built automatically (80% of the residues) is available in

lyso-Gd_SAD/LysoGd_SAD_3_buccaneer.pdb

This model and the experimental map are good enough to start manual building.

Alternatively, run Refmac for a first refinement and manual rebuilt with


LysoGd_SAD_3_buccaneer.pdb

as pdb input file.

Quick SAD phasing in ccp4i2

From the directory `ccp4i2_SAD`, launch `ccp4i2` at the command line
`ccp4i2`


Then **Start a new crystallography project** and define a project directory.



Run the SHELX automated pipeline


In the present section, the automated pipeline for SAD phasing is described. This pipeline is based on SHELX and Refmac5/Buccaneer.

From the `ccp4i2` window, in job list, select **Experimental phasing**, and then the **SHELX** pipeline.



Then provide

- the sequence
 - file **193L.fasta**, to be renamed as .seq
 - directory **lyso-Gd_15May2013/ccp4_MR**
- the crystallography data
 - file **temp_ccp4_f.mtz**
 - directory **lyso-Gd_15May2013/xds_step3_XSCALE-XDSCONV**



and mention that the data were collected at the peak of fluorescence.



Then press **Run** to chain the following steps:

- Substructure determination (SHELXD)
- Phasing, Density Modification, first model (SHELXE)
- Refinement (Refmac5) / Model building (Buccaneer)

Results analysis


Substructure determination:

The substructure consists in 4 Gd atoms, 2 of them with >10 sigma density in the 2mfo-Dfc density map upon refinement, at the interface between 2 molecules, close to residues 62-63.




Model building and refinement:

- 3 cycles of refinement, with final FOM of 0.818, Rfree of 37.5% and Rfactor of 32.3%
- 114 residues built automatically.




Results displayed in ccp4i2 for the model building/refinement step.

The 2mfo-Dfc map, sub-model and refined model can be displayed with coot (press the **Manual COOT** button).



Partial model displayed in the 2mfo-Dfc contoured at 1 sigma. Extra density is visible around the chelating HPDO3A for manual building.



Refined model displayed in the 2mfo-Dfc contoured at 1sigma.

Automated SAD phasing with Phaser

Instead of SHELXE, Phaser can be used for phasing, using the substructure established by SHELXD.

To do so, in the **Experimental phasing** menu, select **SAD phasing from heavy atom sites – Phaser**.

Provide the requested information (wavelength: 1.711888 Ang; MW: 14300).

File
Edit
History
Utilities
Projects
Help

Task menu
View in Coot
View in CCP4mg
Export MTZ
Help
Bibliography
Clone job
Run
Run on server

Job list
Project directory

Job/File	Evaluation
3 SAD phasing from heavy atom sites - PHASER	
1 SHELX	R=0.33 Rfree=0.38

Job 3: SAD phasing from heavy atom sites - PHASER
The job is Pending

Input
Results
Comments

Input data
Keywords

Job title
SAD phasing from heavy atom sites - PHASER

Use data from job
1 Automated structure solution - SHELXC/D/E phasing at
as input below..

Experimental data

Reflections
1 temp_ccp4_f: unknown140515 imported by job 1

Resolution range
Wavelength
1.5418

Known sites, phasing model, or phases

Provide sites, partial model or phases as:
Partial heavy atom substructure

Partial HA model
1 Model coordinates

Composition

For estimating asymmetric unit contents:
Provide an estimate of the molecular weight of protein and

Molecular weight (Da) of protein in the ASU

Molecular weight (Da) of nucleic acid

Cycles of heavy atom model completion
1

Search for pure anomalous scatterers

Show list
Elements to search for

SAD phasing, step by step, with ccp4i

Launch

ccp4i

and define a new project with


mosflm

(the directory in which you have processed the data) as working directory

Step 1: Calculate the anomalous difference Patterson map and check the presence of peaks.

$$P_H(\vec{u}) = \sum_{h,k,l} (|F_{PH}(\vec{s})| - |F_{PH}(-\vec{s})|)^2 \exp[-2i\pi \vec{u} \cdot \vec{s}]$$

Go to menu / sub-menus: Experimental Phasing / Heavy Atom Location / Generate Patterson Map



select “anomalous difference Patterson”

MTZ in: **ctruncate_A5_1_00002-unique.mtz** or **Lyso-Gd_SAD.mtz**

check that F1 stands for F+ and F2 for F-

In “Exclude Reflections” you can choose the resolution range, indicate a threshold in the ratio signal-to-noise, etc.

The list of peaks in the anomalous difference Patterson is stored in **Lyso-Gd_SAD_peaks1.pdb** (Peak coord).

Edit this file with any text editor:

```

CRYST1 77.170 77.170 38.360 90.00 90.00 90.00 P 4/m m m
SCALE1 0.012958 -0.000000 -0.000000 -0.000000
SCALE2 -0.000000 0.012958 -0.000000 0.000000
SCALE3 0.000000 -0.000000 0.026069 -0.000000
ATOM 1 OW WAT X 1 0.000 0.000 0.000 137.19 137.19 O
ATOM 17 OW WAT X 2 25.077 0.000 14.622 9.35 9.35 O
ATOM 13 OW WAT X 3 38.585 35.227 10.790 9.16 9.16 O
ATOM 14 OW WAT X 4 38.585 10.136 12.893 6.86 6.86 O
ATOM 12 OW WAT X 5 15.234 13.308 9.606 6.63 6.63 O
ATOM 12 OW WAT X 6 13.308 15.234 9.606 6.63 6.63 O
ATOM 6 OW WAT X 7 2.715 3.095 4.514 6.40 6.40 O
ATOM 18 OW WAT X 8 25.208 25.208 16.120 6.22 6.22 O
ATOM 4 OW WAT X 9 23.342 23.342 3.180 5.87 5.87 O
ATOM 9 OW WAT X 10 38.585 36.598 6.570 5.72 5.72 O
ATOM 20 OW WAT X 11 28.584 1.994 19.180 5.71 5.71 O
ATOM 7 OW WAT X 12 12.806 9.762 5.186 5.33 5.33 O
ATOM 7 OW WAT X 13 9.762 12.806 5.186 5.33 5.33 O
ATOM 16 OW WAT X 14 11.822 15.680 14.051 5.29 5.29 O
ATOM 16 OW WAT X 15 15.680 11.822 14.051 5.29 5.29 O
ATOM 8 OW WAT X 16 29.094 29.094 5.656 5.27 5.27 O
ATOM 21 OW WAT X 17 21.835 3.112 19.180 5.26 5.26 O
ATOM 3 OW WAT X 18 26.280 26.280 1.214 5.25 5.25 O
ATOM 11 OW WAT X 19 9.264 12.258 9.493 4.64 4.64 O
ATOM 11 OW WAT X 20 12.258 9.264 9.493 4.64 4.64 O
ATOM 5 OW WAT X 21 38.585 16.886 3.847 4.64 4.64 O
ATOM 10 OW WAT X 22 13.464 36.192 8.329 4.62 4.62 O
ATOM 10 OW WAT X 23 36.192 13.464 8.329 4.62 4.62 O
ATOM 19 OW WAT X 24 22.978 28.723 17.844 4.23 4.23 O
ATOM 19 OW WAT X 25 28.723 22.978 17.844 4.23 4.23 O
ATOM 15 OW WAT X 26 26.327 26.327 13.383 3.52 3.52 O
ATOM 2 OW WAT X 27 15.437 38.585 0.000 3.18 3.18 O
END

```

You should have peaks above 4 sigmas.

Step 2: determine the position of Gd atoms using Shelxd program. Basically calculates the coordinates of Gd atoms from the anomalous difference Patterson map.

Go to menu / sub-menus: Experimental Phasing / Heavy Atom Location / ShelxS - Heavy Atom Search

Shelx - Heavy Atom Search Initial parameters from /Espace-no-TimeMachine/Phel...

NOTE This task uses the Shelx program which is not distributed with CCP4 - click Help for details

Job title **ShelxS 2**

Try to find heavy atoms by **Patterson search**

Input format is **MTZ file** using **anomalous** data

MTZ in **Full path..** **/Espace-no-TimeMachine/Phelma-2016/16-11-30-TP-Phelma** **Browse** **View**

FP **FP** Sigma **SIGFP**

DP **DANO** SIGDP **SIGDANO**

Cell Parameters

Space group **P 43 21 2** and lattice type **primitive** Wavelength **1.54178**

Set cell a **77.17** b **77.17** c **38.36** alpha **90.0** beta **90.0** gamma **90.0**

Cell ESDs a **0.001** b **0.001** c **0.001** alpha **0.0** beta **0.0** gamma **0.0**

Search for **4** atoms of **Gd**

Exclude Reflections

Shelx Patterson Search Parameters

Try **10** superposition vectors

Enter 'known' vectors:

Edit list **Add Vector**

Run **Save or Restore** **Close**

MTZ in: **ctruncate_A5_1__00002-unique.mtz** or **Lyso-Gd_SAD.mtz**

check that FP stands for FP and DP for DANO

Guess how many Gd you could have bound to one lysozyme molecule (You can slightly overestimate), and specify that you are looking for Gd.

10 superpositions vectors (number of Patterson peaks that will analysed as Gd-Gd interatomic vectors, you can try from 4 to about 10, not critical in normal cases)

Look at the log file (old style): Two Gd found:

```
Patterson vector superposition minimum function for E ShelxS 2

Patt. sup. on vector 1      0.1222  0.1606  0.2497  Height   6.   Length 18.28

Maximum =      8.96,  minimum =   -98.87      highest memory used = 15501 /673654

      0.1 seconds CPU time

      500 Superposition peaks employed, maximum height  75.1  and minimum height  19.1  on atomic
number scale

Heavy-Atom Location for E ShelxS 2

      9562 reflections used for structure factor sums

Solution   1      CFOM =   1.35      PATFOM = 99.9      Corr. Coeff. = 11.6      SYMFOM = 99.9

Shift to be added to superposition coordinates:  0.4209  0.0618  0.0471

Name  At.No.   x       y       z      s.o.f.  Minimum distances / PATSMF (self first)
GD1   68.4   0.1857  0.5136  0.7047  1.0000  22.47
                                           3.0
```

GD2	63.1	0.1403	0.4813	0.8238	1.0000	18.18	6.27
						2.5	2.7

Maximum memory for Patterson interpretation = 15161 / 4800

0.1 seconds CPU time

Step 2 (alternative)

Use already known phases to locate Gd atoms.

advantages: more precise localization of anomalous scatterers (especially when you have several atoms to locate), no need for very accurate phases

drawbacks: need for phases (either from molecular replacement or from other heavy atoms derivatives)

Here, we can use the phases from the molecular replacement solution done with the bovine lysozyme model (lyso-gd_3.1.mtz contain phases information from the molecular replacement solution)

Go to menu / sub-menus: Map & mask Utilities / Run fft - Create Map

PHIC are the phases from molecular replacements

FOM is the Figure of Merit, i.e. a value between 0 and 1 that quantify the quality of the phase estimates (used here has a weight)

DANO is the anomalous difference $F^+ - F^-$, associated with the sigma (error) SIGANO

The pdb file is only used to determine where the map is calculated (around the model, easier for interpretation, but not mandatory, the map can be calculated on the asymmetric unit if no model is known)

The map calculated is :

$$\rho_H(\vec{r}) = \sum_{h,k,l} w \cdot (|F_{PH}(\vec{s})| - |F_{PH}(-\vec{s})|) \exp[i(\varphi_{\text{calc}} - \frac{\pi}{2})] \cdot \exp[-2i\pi \vec{r} \cdot \vec{s}]$$

with $(|F_{PH}(\vec{s})| - |F_{PH}(-\vec{s})|) = \text{DANO}$; $w = \text{FOM}$; $\varphi_{\text{calc}} = \text{PHIC}$

the shift of $-\pi/2$ come from the fact that we are looking at the f'' contribution (the one that breaks Friedel's law), which contributes with a phase shift of $\pi/2$ (f'' in the structure factor formula)

The highest peaks of this map should correspond to Gd atoms.

The file lyso-1_7_peaks.ha provides the peaks in fractional coordinates:

```
GRID 120 120 64
CELL 77.1700 77.1700 38.3600 90.0000 90.0000 90.0000
ATOM1 Ano -0.1838 0.5129 0.0428 29.37 0.0 BFAC 20.0
ATOM2 Ano -0.3162 0.0129 -0.2928 29.37 0.0 BFAC 20.0
ATOM3 Ano -0.3162 0.0129 0.7072 29.37 0.0 BFAC 20.0
ATOM4 Ano -0.0129 0.3162 -0.2072 29.37 0.0 BFAC 20.0
ATOM5 Ano -0.0129 0.3162 0.7928 29.37 0.0 BFAC 20.0
ATOM6 Ano 0.1838 0.4871 0.5428 29.37 0.0 BFAC 20.0
ATOM7 Ano -0.1403 0.4817 -0.0746 27.04 0.0 BFAC 20.0
ATOM8 Ano 0.1403 0.5183 0.4254 27.04 0.0 BFAC 20.0
ATOM9 Ano 0.0183 0.3597 -0.3246 27.04 0.0 BFAC 20.0
ATOM10 Ano 0.0183 0.3597 0.6754 27.04 0.0 BFAC 20.0
ATOM11 Ano -0.1421 0.4820 0.8906 5.75 0.0 BFAC 20.0
ATOM12 Ano -0.1315 -0.0038 -0.2092 4.02 0.0 BFAC 20.0
ATOM13 Ano -0.1315 -0.0038 0.7908 4.02 0.0 BFAC 20.0
ATOM14 Ano 0.1315 0.0038 0.2908 4.02 0.0 BFAC 20.0
ATOM15 Ano 0.0038 0.1315 -0.2908 4.02 0.0 BFAC 20.0
ATOM16 Ano 0.0038 0.1315 0.7092 4.02 0.0 BFAC 20.0
ATOM17 Ano 0.2835 0.5146 0.2851 3.88 0.0 BFAC 20.0
ATOM18 Ano 0.0146 0.2165 0.5351 3.88 0.0 BFAC 20.0
```

Two independent sites (the first 6 ones are related by crystallographic symmetry and have the same heights (29.37), same for the next 4 ones (heights of 27.04)

the file lyso-gd-fft.pdb provides the peaks in cartesian coordinates:

```
CRYST1 77.170 77.170 38.360 90.00 90.00 90.00 P 43 21 2
SCALE1 0.012958 -0.000000 -0.000000 -0.000000
SCALE2 -0.000000 0.012958 -0.000000 0.000000
SCALE3 0.000000 -0.000000 0.026069 -0.000000
ATOM 5 OW WAT X 1 -14.183 39.579 1.643 29.37 29.37 O
ATOM 5 OW WAT X 2 -24.402 0.994 -11.233 29.37 29.37 O
ATOM 5 OW WAT X 3 -24.402 0.994 27.127 29.37 29.37 O
ATOM 5 OW WAT X 4 -0.994 24.402 -7.947 29.37 29.37 O
ATOM 5 OW WAT X 5 -0.994 24.402 30.413 29.37 29.37 O
ATOM 5 OW WAT X 6 14.183 37.591 20.823 29.37 29.37 O
ATOM 3 OW WAT X 7 -10.831 37.174 -2.860 27.04 27.04 O
ATOM 3 OW WAT X 8 10.831 39.996 16.320 27.04 27.04 O
ATOM 3 OW WAT X 9 1.411 27.754 -12.450 27.04 27.04 O
ATOM 3 OW WAT X 10 1.412 27.754 25.910 27.04 27.04 O
ATOM 23 OW WAT X 11 -10.962 37.199 34.164 5.75 5.75 O
ATOM 6 OW WAT X 12 -10.151 -0.292 -8.025 4.02 4.02 O
ATOM 6 OW WAT X 13 -10.151 -0.292 30.335 4.02 4.02 O
ATOM 6 OW WAT X 14 10.151 0.292 11.155 4.02 4.02 O
ATOM 6 OW WAT X 15 0.292 10.151 -11.155 4.02 4.02 O
ATOM 6 OW WAT X 16 0.292 10.151 27.205 4.02 4.02 O
ATOM 13 OW WAT X 17 21.879 39.714 10.936 3.88 3.88 O
ATOM 13 OW WAT X 18 1.129 16.706 20.526 3.88 3.88 O
```

In order to avoid symmetry related peaks, you can calculate the map on the asymmetric unit:

and the files lyso-1_8_peaks.ha lyso-gd-fft-au.pdb only contains two peaks much higher than the rest (which is basically noise):

```

GRID 120 120 64
CELL 77.1700 77.1700 38.3600 90.0000 90.0000 90.0000
ATOM1 Ano 0.8162 0.5129 0.0428 30.70 0.0 BFAC 20.0
ATOM2 Ano 0.4817 0.8597 0.0746 28.26 0.0 BFAC 20.0
ATOM3 Ano 0.4962 0.6315 0.0408 4.20 0.0 BFAC 20.0
ATOM4 Ano 0.9854 0.7835 0.0351 4.05 0.0 BFAC 20.0
ATOM5 Ano 0.6936 0.3170 0.0308 4.01 0.0 BFAC 20.0
ATOM6 Ano 0.2161 0.3438 0.0380 3.96 0.0 BFAC 20.0

CRYST1 77.170 77.170 38.360 90.00 90.00 90.00 P 43 21 2
SCALE1 0.012958 -0.000000 -0.000000 -0.000000
SCALE2 -0.000000 0.012958 -0.000000 0.000000
SCALE3 0.000000 -0.000000 0.026069 -0.000000
ATOM 16 GD GD X 1 62.987 39.579 1.643 30.70 30.70 O
ATOM 23 GD GD X 2 37.174 66.339 2.860 28.26 28.26 O
ATOM 17 GD GD X 3 38.293 48.736 1.565 4.20 4.20 O
ATOM 11 GD GD X 4 76.041 60.464 1.346 4.05 4.05 O
ATOM 9 GD GD X 5 53.523 24.464 1.182 4.01 4.01 O
ATOM 10 GD GD X 6 16.678 26.533 1.457 3.96 3.96 O

```

Here, the peaks of the map gives you directly the coordinates of the Gd atoms.

Further analysis should show that these two peaks match the two peaks found by difference Patterson map, if you take into account the possible symmetries.

Step 3: calculating the experimental phases using the anomalous signal and the two gadolinium sites determined previously.

A) the lyso-gd-fft-au.pdb needs some editing to be interpreted properly (it will be saved in lyso-gd-fft-aus.pdb):

- just keep the two main sites
- put the occupancy to 1.0 (it will be refined anyway afterwards)
- replace o by GD in the last column (otherwise the program will believe you have oxygen instead of Gd, and will fail)

```
CRYST1    77.170    77.170    38.360   90.00   90.00   90.00 P 43 21 2
SCALE1      0.012958 -0.000000 -0.000000      -0.000000
SCALE2     -0.000000  0.012958 -0.000000      0.000000
SCALE3      0.000000 -0.000000  0.026069     -0.000000
ATOM       16  GD   GD X    1      62.987  39.579   1.643   1.00 30.70      GD
ATOM       23  GD   GD X    2      37.174  66.339   2.860   1.00 28.26      GD
END
```

B)

Go to menu / sub-menus: Experimental Phasing / Automated Search & Phasing / Phaser SAD pipeline:


Pay attention to deselect run SHELXD before Phaser (you have already found where the GD atoms are)

some files are in the ccp4 directory and not in the mosflm one: just be sure where they are in you session.

Here you need F+ and F- and not DANO

The main 3 steps of the program will be:

- 1) Using Phaser to calculate phases from the data and positions of Gd atoms, using the principles described in the lectures. Positions, B-factor and occupancy of Gd atoms will be refined. Other sites will be searched and refined. (check log file)
- 2) Using Parrot for density modification procedure (solvent flattening, etc ...): should improve the quality of the phases
- 3) Using Buccaneer to build automatically a model in the resulting electron density map, using the primary sequence information provided.



the experimental phases are in lyso-1_10_parrot.mtz. The model automatically constructed by buccaneer is in lyso-1_10_buccaneer.pdb

Files with 'hand' are also created in cases where the correct hand is not known. Here, the space group was set properly, since it was already known.

Step 4: Use coot to look at maps and models

open coot

open the molecular replacement solution with bovine lysozyme in ccp4 directory: lyso-gd_3.1.pdb
menu File / Open coordinates

open molecular replacement maps in ccp4 directory: lyso-gd_3.1.mtz
menu File / Auto Open MTZ

open the automatically built model in mosflm directory: lyso-1_10_buccaneer.pdb
menu File / Open coordinates

Open experimental phases map in mosflm directory: lyso-1_10_parrot.mtz
menu File / Open MTZ, mmCIF, fcf or phs... and chose
parrot.F_phi.F as amplitudes
parrot.F_phi.phi as phases

adjust the contour to 1 sigma

Compare the different maps. Focus on regions where the sequence of bovine and Hen lysozyme are different.

Hen and bovine lysozymes sequence alignment :


```
>2z2f_A Lysozyme C-2; stomach lysozyme, 1,4-beta-N-acetylmuramidase C, bacteriolytic enzyme, hydrolase; 1.50A {Bos taurus} SCOP: d.2.1.2
```

```
Length = 129      Score = 154 bits (389), Expect = 1e-38,      Method: Composition-based stats.  
Identities = 71/130 (54%), Positives = 96/130 (73%), Gaps = 2/130 (1%)
```

```
Query: 19  KVFGRCELAAAMKRHGLDNYRGYSLGNWVCAAKFESNFNTQATNRNTDG-STDYGILQIN 77  
          KVF RCELA  +K+ GLD Y+G SL NW+C  K+ES++NT+ATN N   STDYGI QIN  
Sbjct: 1   KVFERCELARTLKKLGLDGYKGVSLANWLCLTKWESSYNTKATNYPSSSESTDYGIFQIN 60
```

```
Query: 78  SRWWCNDGRTPGSRNLCNIPCSALLSSDITASVNC AKKIVSDGNGMNAWVAWRNRCKGTD 137  
          S+WWCNDG+TP + + C++ CS L+ +DI +V CAK IVS+  G+ AWWAW++ C+  D  
Sbjct: 61  SKWWCNDGKTPNAVDGCHVSCSELMENDIAKAVACAKHIVSE-QGITAWVAWKSHCRDHD 119
```

```
Query: 138 VQAWIRGCRL 147  
          V +++ GC L  
Sbjct: 120 VSSYVEGCTL 129
```



MR

under construction

Complete model

Partial/remote model

Refinement / construction

under construction

Complete model

Partial/remote model